


User Manual

Bio Matrix Scan Software



Please read the manual before using the tools.

Version 4.2

Contents

1. Introduction.....	3
Features.....	3
2. Package.....	4
3. Description.....	5
Warranty.....	5
4. Procedures.....	6
Software Installation.....	6
The Program Interface.....	8
Survey Result Tabs (Brief Details)	9
5. The Program Operation.....	11
Program Start.....	11
Adding A New Client.....	12
Removing A Client from The List.....	13
ECG and ECG Tab.....	13
Rhythmogram.....	14
Dynamics of Physiological Parameters During The Survey.....	14
Functional status Indicators.....	15
6. Survey Result Tab (In Details)	16
6.1 The Vegetative Regulation Tab.....	16
Histogram of «R-R» intervals.....	17
6.2 The Neurohumoral Regulation Tab.....	18
6.3 The Psychoemotional State Tab.....	20
6.4 The Fractal Analysis Tab.....	22
6.5 The Complex Analysis Tab.....	26

6.6 The Dynamics of Functional Condition Tab.....27

6.7 Portrait of The Aura.....29

6.8 Chakras Activity Map.....30

6.9 The Chart of The Meridians.....33

7. Save The Report.....39

8. Printout of Survey Reports.....40

9. Change Program Language.....41

10. Effectiveness of Treatment Evaluation.....42

Long-term Observation.....42



1. Introduction

Bio Matrix Scan (BMS) is intended to be used for both assessing the state of one's health and controlling the efficacy of treatment. BMS is distinct from other systems in that it is simple to use. Operating the unit does not necessitate any particular medical training. It is possible to obtain survey results fast and easily by following the complete step-by-step user instructions.

BMS allows for the visual depiction of the body's reaction to medications and medical procedures.

Features

- Assess the status of the cardiovascular and vegetative nerve systems using variance analysis.
- Determine the degree of organism adaption and biological rhythm harmonization using the fractal analysis method.
- The user's biological age is determined.
- Displaying the results of two surveys at the same time in order to do a comparative analysis.
- The survey's results can be printed.

2. Package



- | | |
|--------------------------|----|
| 1. Bag | x1 |
| 2. Right-sided Electrode | x1 |
| 3. Left-sided Electrode | x1 |
| 4. Bio Matrix Scan | x1 |
| 5. USB-C | x1 |

3. Description

Warranty

Each single product passes through a quality process. Nevertheless, if failures occur, contact the customer service at once. The product is under 12 months (1 year) warranty effectively according to the delivery date stated on the warranty card.

Warranty does not apply to damages, which result from misuse or neglect. The warranty also expires if you dismantle the product or failed to provide the warranty card. The manufacturer is not liable for consequential damage or in case of a non-intended use of the product.

If a failure occurs during the warranty period, the product will be replaced or repaired. The manufacturer reserves the right to exchange components of the product instead of repairing it. If the failure results from misuse or neglect, the user must pay for all cost that may occur. The freight cost and other costs will be depending on which Term & Conditions criteria the product falls into. **Import tax and Customs clearance duty cost is not included in the term of warranty.** In that case, you may ask for an estimated cost beforehand.

4. Procedures

Software Installation

Download Software from this link:

<https://lmap.de/emailclients/MandalaSoftwareGen2.2021-10-04.zip>


- Double click the zip file.




- When accessing to the zip file, double click the file.




- Select a language and then do the steps accordingly.





■ Click to finish installing

The Program Interface



Main Interface

Survey Result Tabs (Brief Details)


Displaying a survey on a panel. The findings for the current client are displayed on the right-hand side of the window. This panel consists of a series of tabs. The buttons at the top of the screen are used to switch between them.



There are eleven tabs:




- The control elements for recording and displaying the client's electrocardiogram are found on the Registration and Viewing of Data page.



- The Vegetative Regulation tab shows how heart rhythm variance analysis is used to determine the organism's vegetative regulation.



- The Neurohumoral Control tab displays a neurodynamic analysis-based hormonal regulation evaluation.



- The Psychoemotional State tab shows the measurement of psychoemotional state using brain biorhythm mapping.



- The Fractal Analysis tab shows the body's adaptation level and is used to determine biological age by fractal analysis.



- The Complex Analysis tab displays a complex analysis of the organism's functional state.



- The Dynamics of Functional Condition tab shows how the change in functional state indices changes over time from results.



- Portrait of the Aura is an examination of energy processes within the human body, processes that are expressed in a heart rate that objectively reflects the cardiovascular, nervous, endocrine, and vegetative systems' quality.



- The aim of the Chakras Activity Map is to examine energy points in your body that correspond to nerve bundles, major organs, and areas of our energetic body that influence our emotional and physical well-being.



-The Chart of The Meridians is the chart of the analysis of the twelve major meridians.



-Light Mandalas' Copyright Information.


5. THE PROGRAM OPERATION


Program Start

After software installation, to begin, you can plug in BMS to one of your Laptop's USB sockets.

Then, move the mouse pointer over "Mandalas" icon  on the desktop and double-click with the left mouse button.

The application window will appear on the screen once the program has loaded:



The green lamp in the bottom right corner of the primary window must be turned on , indicating that the unit is ready to use.





In brackets number of virtual USB port is specified. Color of a bulb means readiness of the device for work:


- **Gray** – The device is not found.
- **Red** – It is wrong connected device.
- **Green** – The device is ready to work.

The button  allows turning on or off sound signals of the program.

Adding a new client

A screenshot of a software interface for adding a new client. The interface is a vertical form with a light gray background. At the top is a large white rectangular area for a profile picture, containing a gray silhouette of a person. Below this are three input fields: the first is labeled 'Patient-1' and has a blue male sex icon to its right; the second is labeled '2/15/1996' and has a calendar icon to its right; the third is labeled 'E-Mail'. At the bottom of the form is a small, faint icon.

Click icon "Add" in the bottom left of the main interface,  Add. A new client will be added to the List of clients with an empty personal information area. Enter the added person's first and last name in the "Name" area. The client sex must also be entered. To do so, click on the  when you press the male sex icon with your left cursor, it will change to a female sex icon  after each iteration of clicks, then back to the male sex icon. .

The date of birth of the newly added client should be entered in the field "Birth Date." The date can be entered numerically (in "DD/M/YY" format) or by using the calendar. .

To add the client's photo, press the "the client's Photo" area with the left mouse button. The regular Windows panel will open, allowing you to pick the client's photo file. Any image file, regardless of its size, can be used as the client's picture.

Removing the client from the List


Pick the client that should be removed from the List of Clients and click “Delete”

Delete

with the left mouse button. The window of client removal confirmation will appear if the selected client has any survey data. If no surveys have been completed, the client will be deleted from the system without any more questions of confirmation.

ECG TAB


Brief description of control elements

-  /  - ECG recording start/stop button.
-  button of polarity change of recorded ECG signal.

 89 - heart rate.

R-R 300 - number of recorded heart contractions.

ECG



The software shows either an ECG signal recorded previously during the survey or an ECG signal recorded immediately at the moment, depending on the selected operation mode. The ECG amplitude is laid off along the vertical axis, while the time in minutes and seconds elapsed from the registration start is laid off along the horizontal axis. The graph scaling is done by rotating the mouse while pressing the right button. Dragging the graph is achieved by dragging the mouse while holding the left button.

Rhythmogram



The software shows either a previously recorded rhythmogram or the one currently being registered, depending on the activity mode selected. The number of RR intervals is plotted on the horizontal axis, while the interval duration in seconds is plotted on the vertical axis in a rhythmogram. The RR interval is the amount of time between two heartbeats. Objects, such as extra systoles and disruptions, are outlined in red on the rhythmogram. The right mouse button is used for graph scaling, while the left mouse button is used for graph dragging.

The button  allows you to see the Graph of Functional Condition Dynamics During Recording.

Dynamics of Physiological Parameters During the Survey

The Graph of the Dynamics of Physiological Parameters During the Survey shows how the functional state indices changed during the survey. It allows the survey results' reliability to be estimated.

If the usable state remains constant during the study, it indicates that the survey was properly performed, and the findings are trustworthy. If the graph of functional state, on the other hand, has sudden level decreases, it suggests that either the client was not at rest during the survey, or that the ECG signal received from the client was influenced by external disturbances.



If the results of such a survey are unreliable, so a repeat trial should be conducted after removing the variables that have a negative effect on the consistency of the registered ECG signal.

Functional Status Indicators



The indicators display the normalized values of functional status indices:

“A” – level of organism adaptation

“B” – vegetative regulation index

“C” – neurohumoral regulation index


“D” – psycho-emotional state index

“H” – Health, complex state index

6. Survey Result Tabs (In Details)

6.1 The Vegetative Regulation Tab

Vegetative regulation which controls physiological processes below the level of human consciousness, performs vegetative control. It responds quickly to changes in the internal and external environment by affecting the cardiovascular system, whose proper functioning controls the supply of nutrients and oxygen.






Histogram of «R-R» intervals

It is a graph that depicts the relationship between the number of «R-R» intervals and their length at various time intervals. The length of «R-R» intervals is plotted along the horizontal axis, while the likelihood of their occurrence (i.e., the number of «R-R» intervals within the respective range) is plotted along the vertical axis. The interval between histograms is 0.04 seconds. The histogram base ranges from 0.32 to 1.64 seconds. The central location of the bars, with the highest bar localized within the range of 0,7–1,0 s, characterizes vegetative equilibrium. A major change to the left, as well as a narrowing of the histogram base, is seen when the sympathetic division of the vegetative nervous system is dominant. When parasympathetic influence is present, the opposite effect is seen.

Tension index: The degree of myo-card tension is defined by this variable.




Autocorrelation portrait expresses the degree to which different rhythm gram fragments are identical.

Vegetative balance index establishes a connection between the sympathetic and parasympathetic divisions of the vegetative nervous system's function.


6.2 The Neurohumoral Regulation Tab

The Neurohumoral Regulation Tab shows data about the neurohumoral control parameters.


Neurohumoral regulation system maintains the consistency of the internal environment and the organism's long-term adaptation to evolving living conditions by regulating the composition and structure of biochemical substances in the organism.



The structure of physiological rhythms is described as a neurodynamic matrix, with each variable describing the dynamics of the corresponding rhythm. Individual matrix elements represent the rhythms of the corresponding organism systems, with the color of each element indicating how closely the parameters of these rhythms match the universal law of living nature behavior, namely the law of two exponents. The parameters of the "ideal exponent" are in line with the "golden ratio." Compliance with these parameters ensures that life-support systems operate as well as possible while using the least number of resources possible. The yellow and red colors of matrix components show that the rhythm's parameters are not optimum.




Neurohumoral Regulation Index identifies the efficiency with which the nervous system functions and how efficiently the organism utilizes its energy and physiological resources. The neurohumoral control mechanism is in charge of maintaining the internal environment's consistency and adapting the organism to changing living conditions.



Energy Pyramid describes the total amount of an organism's physiological resources as well as the equilibrium between resource accumulation and consumption cycles under current life speed conditions. The dynamics of anabolic and catabolic processes in the body are characterized by the correlation between the left and right parts of the "pyramid." The blue area of the pyramid represents the time it takes to recover resources, while the red area represents the time it takes to consume them. The smallest scale of the entire pyramid indicates that physiological resources have been depleted.

6.3 The Psycho-emotional State Tab

Psychoemotional State Tab displays information on the parameters of client's psychoemotional state.




The **spline-map** is the product of spline-interpolation of the main brain rhythms' peak values. Spline maps of brain electrical activity are a series of bivariate distributions of central nervous system (CNS) rhythms in the brain's functional spaces that characterize the CNS's integral activity in these spaces.



The hypoactivity of the brain as a result of illness, overfatigue, and stress is shown by areas of black and yellow-red color.

Psychoemotional state index describes the extent to which stress has a negative impact on the body.

The frequency **spectrum of brain activity** is obtained by translating heart rhythms into brain rhythm frequency ranges. The obtained spectrum is divided into frequency ranges for delta, theta, alpha, beta, and gamma rhythms.



The normal state is characterized by an even distribution of rhythms across the entire frequency spectrum.

Delta rhythm predominance is suggestive of illness, exhaustion, and stress.

Theta rhythm predominance is suggestive of meditation state, calmness and concentration.

Alpha rhythm predominance is suggestive of relaxation.


Beta rhythm predominance is suggestive of logical thinking state and decision-making.

And Gamma rhythm predominance is suggestive of future telling and higher power spirituality.

6.4 The Fractal Analysis Tab

Details regarding the **Fractal Analysis** of the client's status is shown in the Fractal Analysis Tab


Fractal analysis is used to visually determine the degree of harmonization of biorhythms in fractal-structured body organs and structures in order to identify functional and pathological changes, assess the immune status of the organism, and predict client health condition changes.



Fractal analysis reveals regularities in complex parameters of rhythms that describe physiological processes in the human body. Rhythms with a fractal structure that manifest as universal reactions to external and internal stimuli are identified during electrocardiogram analysis. The degree to which these rhythms are in sync defines the quality-of-life behavior at any given time.

The ability of an organism to sustain and retain this balance over time reflects its ability to respond to evolving external and internal environmental factors, or its "immunity."

The fractal portrait index is a data index of an organism's immune status that demonstrates its ability to adapt to new situations.



Mandelbrot's equation, whose parameters are the dynamic parameters of ECG, is represented by a fractal portrait. The degree of harmonization of physiological rhythms is shown by a fractal portrait, which shows how well the organs and systems of the body work together.

Fractal harmony of biorhythms at all levels of control suggests abundant energy supplies, an optimum energy supply balance, and a favorable prognosis for health state improvements.


Partially disrupted biorhythms in some organs and systems of the body may indicate functional disturbances in those organs and systems.

Disturbance of biorhythms' fractal equilibrium, which leads to severe functional disruptions or pathological changes in individual organs, characterizes the depletion of energy resources and degradation of energy balance.

Level of harmonization of biorhythms characterizes the consistency of regulatory systems, represents the efficiency of an organism's adaptability to evolving circumstances, and serves as a knowledge index for immune system health.

The **gerontological curve** represents the biological rhythm that is the slowest. It has a lifespan of 90-110 years, which is determined by life expectancy. It's a graph that combines logarithmic and exponential curves to show the relationship between the rates of life resource accumulation and consumption.

These correlations correspond to the parameters of the "golden segment" in the case of the "model" gerontological curve. The nominal "phase" at which the accumulation of life resources transitions to their slow consumption is estimated to be about the age of 30-33 calendar years.



The clock on the outside of the body is used to calculate calendar age, while the clock inside the body is used to calculate biological age.

People may be younger or older than their chronological age, and as a result, illnesses and deaths associated with aging occur at various calendar ages. The majority of these disorders have hereditary causes. As a result, a person's true age is determined not by the passage of time after birth, but by the factors that reflect his or her viability. The biological age of a person refers to the body's viability level. Biological age refers to the systemic disintegration of an organism as it ages.

We may assess a person's biological and true age using the principle of "gerontological curve," which does not always correspond to their calendar age, for better or worse.

The model gerontological curve was created using statistical analysis of biorhythms in over 10,000 clients of various ages. It is a graphic depiction of the accumulation and consumption rates of life resources in the body of an average individual over the course of a 100-year life cycle.

The amount of life resources is proportional to the level of integral health index. We get the point of client's "biological age" by marking this index on the ordinate axis and projecting it on the model gerontological curve. The biological era is precisely the projection of this point on the abscissa axis in the traditional sense of calendar calculation.



6.5 The Complex Analysis Tab

The Complex Analysis Tab displays the information on the complex analysis of client's status.




The Complex State Index is a symbolic statistical representation of a "client's health condition" that combines all the previous indices. It aids in determining the overall consistency and equilibrium of the body's physiological processes. The term "health level" refers to the amount of stress that the body's regulatory mechanisms need to achieve balance with the surrounding world.


The normal for an average client of this age is indicated by the thin blue disc in the Diagram of Complex State Index, while the literally calculated health index can rank above or below the standard level. The contrast between the individual and average value is seen in the diagram's top left corner.

6.6 The Dynamics of Functional Condition Tab

The process of change of functional status indices over time is shown in the Dynamics of Functional Condition Tab.



Two date management elements are located at the top of the tab and can be used to show the analyzed era. By default, only the ten most recent surveys are examined.



The graph "Dynamics of Functional Condition" shows the outcomes of all the examined surveys.

The buttons **A B C D H**  in the graph's upper right corner allow you to hide or show the corresponding parameter.

One of the surveys on the graph can be selected with the left mouse button to see more specific information on it. With a yellow vertical bar, the currently selected survey is highlighted.


The information on the currently selected survey is shown below the graph.



The date and time of the survey, as well as the heart rate during the survey and the values of the functional state indices, are all shown. In addition, a regular prediction of the client's condition is made based on the results of the considered survey study. The accuracy of prediction is directly proportional to the number of surveys examined. Surveys performed more than a month ago are not considered when calculating forecasts.


The calendar in which survey dates are highlighted is located at the bottom of the Dynamics tab. The highlight color represents the client's current condition. If multiple surveys were performed on the same day, the highlight color would equate to the averaged value of the client's condition. By left clicking on the highlighted date in the graph "Dynamics of Functional Condition," you can choose the appropriate survey.

6.7 Portrait of The Aura



"Aura" is nothing more than the ability to sustain and retain the equilibrium of physiological processes over a long period of time, according to science. It reflects the organism's ability to adapt to changing internal and external environments.

In reality, "aura" is a data indicator of the body's immune status, showing its capacity to adapt to new circumstances. The study of energy processes within the human body, such as those expressed in heart rate, enables us to create an integral index that objectively reflects the quality of the cardiovascular, nervous, endocrine, and vegetative systems. From the standpoint of physical substance, this index is like the idea of "aura" as known by oriental healers.

The button  allows for the creation and storing of sound files on a computer disk, which can be used to affect the client's aura and thereby enhance both her condition and the client's overall condition.

The dynamics of the integrated index

6.8 Chakras Activity Map



Description of chakras

1. The Crown Chakra - its colour is Violet/Purple and it's at the top of our head. The Crown Chakra is our sense of Beauty, Creativity & Inspiration. It relates to self knowledge & spiritual awareness and is the union with our higher self. Disease (Dis-ease) can result with an in-balance of energy in this Chakra, either too much, or too little. The purple energy of the Crown connects us to our spiritual self bringing guidance, wisdom & inner strength. It purifies our thoughts & feelings and gives us inspiration in all of our actions. It enhances artistic gifts & creativity.

Associated problems with the Crown Chakra are: Depression, Parkinson's disease, Schizophrenia, Epilepsy, Senile Dementia, Alzheimer's, many mental disorders, confusion and dizziness.

Personality Traits are: Inspirational leadership, kind & fair, humanitarians, visionary, self sacrificing, creative and mentally strong.

2. The Third Eye Chakra - its colour is Indigo and is in the centre of your the forehead. The Third Eye Chakra is our Intuition, Mysticism & Understanding. It is related to the eyes, lower head and sinuses. Indigo relates to responsibility of the self, being responsible for your own life and to follow your soul's path. To trust your own intuition. Indigo energy connects us to our unconscious self (not the thinking self/ego) and give us the experience of being part of the universe. the Third Eye Chakra strengthens imagination, psychic powers & increases our dream activity.

Associated problems with the Third Eye Chakra are: Tension, headaches, migraine, visual defects, short and long-sightedness, cataracts, glaucoma, sinus problems and ear ache/problems.

Personality Traits are: Intuitive, fearless, idealistic, practical, wise and a truth seeker.

3. The Throat Chakra - its colour is blue and it of course located in our throat. The Throat Chakra relates to our Knowledge, Health & Decisiveness. Blue is the colour of spirit and relates to self expression, speech, communication, the ability to communicate our needs, it's the spirit of Truth & Purpose. Blue has a pacifying effect on the nervous system & brings beautiful relaxation which is ideal for insomnia, and hyper-active children. It brings clarity to our speech and gives us wisdom.

Associated problems with the Throat Chakra are: Thyroid problems either over or under active, anorexia, asthma, bronchitis, hearing problems, tinnitus, mouth ulcers, sore throats and tonsillitis.

Personality Traits are: Loyal, tactful, affectionate, inspiring, inventive, caring & cautious.

4. The Heart Chakra - its colour is green and it governs self control, love and balance. It relates to love/self love and the ability to give and take it unconditionally. When we are balanced we are able to give love as well as love and nurture ourselves. Green helps relax our muscles, nerves and thoughts. The heart cleanses & balances our energy to give a feeling of rebirth, peace & harmony.

Associated problems with the Heart Chakra are: Heart disease, AIDS, Chronic Fatigue Syndrome, allergies and breast cancer. Pink also relates to the heart being the colour of love.

Personality Traits are: Understanding, adaptability, sympathetic, compassionate, generous, humane, a lover of nature, romantic and self controlled.

5. The Solar Plexus Chakra - its colour is yellow and is situated just below the ribs. It governs our wisdom, clarity & self esteem. Yellow is a creative colour and it relates to our self worth, how we feel about ourselves and how others view us. This is the area of our personality, ego & intellect. It gives us clarity of thought, increases our awareness and stimulates interest and curiosity. Yellow energy is related to the ability to perceive & understand, it connects us to our mental self.

Associated problems with the Solar Plexus Chakra are: Diabetes, pancreatitis, liver disease, peptic ulcers and gall stones.

Personality Traits are: Good sense of humour, optimism, wisdom, confidence, practical and intellectual.

6. The Sacral Chakra - its colour is orange and it is situated in the lower abdomen. It governs our happiness, confidence and resourcefulness. Orange is the colour of success and relates to self respect; having the ability to give ourselves the freedom to be ourselves and helps us to expand our interests and activities. Orange brings joy to our workday & strengthens our appetite for life. It's the best emotional stimulant! It connects us to our senses and helps to remove inhibitions, makes us social and independent.


Associated problems with the Sacral Chakra are: pre-menstrual syndrome, problems with menstrual flow, uterine fibroids, ovarian cysts, irritable bowel syndrome, endometriosis, testicular disease and prostate disease.

Personality Traits are: Enthusiastic, happy, sociable, energetic, sporty, self assurance and constructive.

7. The Base Chakra - its colour is red and it is situated at the base of the spine. It governs our vitality, courage & self confidence. Related organs are the kidneys, bladder, hips and legs. Red relates to self awareness and is the area of survival and your place on this earth. Red provides the power from the earth and gives us energy on all levels. It connects us to our physical body...everything that is to be commenced needs the vitality of red!

Associated problems with the Base Chakra are: Constipation, diarrhea, piles, colitis, Crohn's disease, cold fingers and toes, frequent urination, high blood pressure, kidney stones, impotence, hip problems, legs and feet. Pink can also be used in this area for a more gentle energy if required.

Personality Traits are: Courageous, confident, humanistic, strong willed, spontaneous, honest and extroverted.



6.9 The Chart of The Meridians

The details about the state of the client's meridians are reflected in the insert "MERIDIAN." The diagram of meridians is shown on the projector. The shape and color of this is determined by the state of the client's meridians. The capacity of each meridian to embrace and give energy during a regular period of meridians is also shown on the computer.



On the screen, relative sizes of movement of the client's meridians are also shown.



Description of human meridians

It is identified as the most inactive meridian in the study of meridional diagnostic data for forecasting of an anticipated target organ. Then his Yin – Yang strength is calculated. It allows you to select competently individual correction methods and check the predicted result in real time.

The general portrait (color, aura intensity, gaps and protuberances), percent of the general, is estimated by visual analysis of the AURO-gram. In a choro-gram study, their behaviors are compared in terms of radiation intensity or a percentage ratio

1. GB – Meridian of the Gallbladder (11P.M.-1A.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	feeling of completeness in a stomach, weight and pains in the head, a pain and a spasm in cubcostal area	disorders of Gall Bladder functions, pains along the Meridian, an impairment of a gastroenteric tract, neuropsychic disorders	to stop injections and drug intake, to address to the attending physician, to refuse alcohol, meat, to avoid physical activities and weariness of eyes
	Lack of energy	breakdown, tearfulness, feeling cold, dizziness, hindered breath, heavy and deep sighs, a lurch gait		

2. LR – Meridian of Liver (1A.M.-3A.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	tearfulness, irritability, irascibility, morbidity, impulsiveness, an aspiration to work without interruption while all work won't be finished	Pains in the area below ribs at right, jaundice, pains along the Meridian and in the waist, diseases of genitals, a cold, fever	to stop injections and drug intake, to address to the attending physician, to refuse alcohol, meat, to avoid physical activities and weariness of eyes
	Lack of energy	a lurch gait, awkwardness, dizziness, fast fatigue		

3. LU – Pair Meridian of Lung (3A.M.-5A.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	Pneumonia, small capacity of lungs, a strong cough with a pain, plentiful phlegm	Diseases of lungs and upper airways with the stagnation phenomena, and also diseases of cardiovascular system	To refuse smoking, to avoid some excessive cold and hot air, to breathe fresh air
	Lack of energy	Algor, tussiculation, feeling cold		

4. LI – Meridian of the Large Intestine (5A.M.-7A.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	Rumbling in a digestive tract, dryness of lips and in a mouth, a lock	restriction of movement and pain in joints of the top extremities, inflammatory processes along the Meridian, frustration of digestive organs, especially the intestines	To refuse smoking, to avoid some excessive cold and hot air, to breathe fresh air
	Lack of energy	A pain and immovability in shoulders and hands, dizziness, diarrhea, sometimes a lock		

5. ST – Meridian of the Stomach (7A.M.-9A.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	overeating and a heat in a body, a pain and spasms on external lateral surface of legs	digestive disorders, a pain along the Meridian	not to allow an overeating, not to take alcohol and sweets, to avoid humidity
	Lack of energy	appetite loss, sometimes a complete refusal of food, a weakness in feet and their cold snap		

6. SP – Meridian of the Spleen and Pancreas (9A.M.-11A.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	unstable appetite, feeling of weight and numbness in a	pains in the area epigastria, below the ribs and in the	not to allow an overeating, not to take alcohol and

		body, an aspiration frequent to have a rest, and to lie down	chest	sweets, to avoid humidity
	Lack of energy	a passion to sweets, memory impairment and drowsiness at the afternoon, meteorism		

7. HT – Meridian of the Heart (11A.M.-1P.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	a heartache, a pain in shoulders and forearms, excitability, feeling weight in breast, sometimes temperature increase and dryness in a mouth	a pain in the heart area, cardiac rhythm disorders, an anxiety, decrease in memory	to increase frequency of food intake, to reduce its quantity, in one reception the food should be same
	Lack of energy	indecision, depression, sensation of fear, excitement, tachycardia, short wind, dizziness		

8. SI – Meridian of the Small Intestine (1P.M.-3P.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	a pain round a neck or in one half of the head, noise in ears, decrease in hearing	disorders of passableness through the intestines, pains along the Meridian, psychological frustration	to increase frequency of food intake, to reduce its quantity, in one reception the food should be same
	Lack of energy	a pain in temples, lateral surface of a neck, forearms, a cold snap of extremities		

9. BL – Bladder Meridian (3P.M.-5P.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	pains in a waist and lower extremities, frequent painful urination with small amount of urine, an eye pain, a headache, spasms in gastrocnemius muscles	Cystitis, cystalgia, nephritis, heart troubles, diseases of eyes, pains in a backbone and a coxofemoral joint, a headache, bleeding from a nose, hemorrhoids	to abstain from sexual life, to avoid physical and intellectual overloads, humidity, excessive heat and cold, to reduce liquid and salt consumption
	Lack of energy	rare plentiful urination, backbone pains, feeling weight in feet and a back, dizziness, decrease in sight		

10. KD – Meridian of the Kidney (5P.M.-7P.M.)

Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	an unusual inflow of forces, an aspiration to work without interruption, noise in ears, dark brown color of urine	urinary disorders, diseases of genitals, neuropsychic disorders, weakening of vital forces and a predisposition to fatigue	to abstain from sexual life, to avoid physical and intellectual overloads, humidity, excessive heat and cold, to reduce liquid and salt consumption
	Lack of energy	decrease in will power and sexual inclination, impatience, uncertainty, a cold feet, intense neck		

11. PC – Pericardium Meridian (7P.M.-9P.M.)


Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	a pain, feeling weight and heat in the head, a belly-	diseases of heart, lungs, mental disorders, pains	habitually used quantity of foods is to lower to 1/3 and to

		ache, a lock, heat in palms	along the meridian	avoid physical and intellectual overloads
	Lack of energy	bad dream with a large number of dreams, short wind, a fear of height, fatigue, desire to lie down, diarrhea		

12. TE – Meridian of Triple Heater (9P.M.-11P.M.)


Power Activity	Power Balance	Signs	Diseases	Recommendations
Less than 30%	excess of energy	weakening of hearing, a pain in the neck, shoulders, shovels, top extremities, ears, bad tolerance of heat, irritability, sleeplessness	arthritises of hands joints, noise in ears, deafness, depressions	habitually used quantity of foods is to lower to 1/3 and to avoid physical and intellectual overloads
	Lack of energy	a temple pain, algor, weakness		

7. Save the Report in PDF file




Click  and then select “print report”

Select “Microsoft Print to PDF”
and then click “OK”




Name PDF file and click “OK”




Scan Result Report Example

8. Printout of Survey Reports


To print the results of any client's survey, use the left mouse button to select the control factor Select Survey Date, then select Print Report from the menu that appears.



Click  and then select "print report"

Select your printer and then click "OK"

9. Change program language



Click 

Select “Change program language”,
then click the language you want.

10.EFFECTIVENESS OF TREATMENT EVALUATION

Initial indices are measured in the client, during which the client is subjected to some type of clinical intervention (fast-acting drug intake, physiotherapy, reflex therapy, etc.) and indices are measured again immediately after the therapy exposure.

The results of repeated surveys are compared to the results of the first survey.

The therapeutic intervention has had a beneficial impact if the functional condition has changed.

The therapeutic effect will initiate the organism's healing process; but, in this situation, indices may deteriorate as the body requires resources to fight the disease. If the indices have declined after the corrective procedure, it is advised to wait for a short period of time (15-30 minutes) before repeating the measurement.

Long-Term Observation

Long-term observation is used in order to keep track of the client's functional state on a daily basis. Before beginning the course, it is important to assess the initial indices, and then perform periodic surveys. It is preferable to perform all surveys at the same time and at frequent intervals in order to collect objective data.

Disclaimer

All content found on the Light-mandalas.com Website, Facebook Page, YouTube, including: text, images, video, audio, or other formats are intended for informational, educational and research purposes only and neither Light Mandalas Co. Ltd. nor the publishers of this content take responsibility for possible health consequences of any person or persons reading or following the information in this educational content and is not intended to be a substitute for professional medical advice, diagnosis, treatment or cure. Always seek the advice of your physician or other qualified health professionals with any questions you may have regarding a medical condition.

“No claims other than stress detection and stress reduction may be made”

Light Mandalas products are to be used as a Heart Rate Variability Biofeedback system and Light, Color and Sound therapy, designed for stress detection and stress reduction only. The devices do not evaluate any issue other than stress. Stress can come from many sources; this system uses light and color and sound therapies to reduce stress. Remember, there has never been any reported side effects using light, color and sound therapy anywhere in the world and throughout time.

Notes on the practices of low-level laser therapy (LLLT) and Biofeedback

Low-level laser therapy (LLLT) is considered investigational by the FDA, and no medical claims of cures are permitted in the US; at this time however, low-level lasers are defined within the scope of complimentary medical practices such as physiotherapy, aesthetic, Traditional Chinese Medicine and acupuncture, etc.

The FDA does not regulate Bio resonance tools but they do regulate biofeedback tools and some Bio resonance tools are registered under this category.

*It is a legal requirement at this time that the user does not make any medical claims of diagnosis or cure.

Precaution

Light Mandalas products are working on energetic level. The standard benefit is energy clearance. During the process, your body and mind can be uncomfortable. The symptoms associated with overly rapid clearance of toxins. If those detox symptoms are too strong, you may need to stop using the products. You have to allow the clearance or detox to proceed and you can start to try the products afterwards. However, if the issues still remain the same, it may need to be completely discontinued for your own safety

Shall there be any inquiry, please contact us for further support.

Light Mandalas Co., Ltd. (Thailand Headquarters):

- **Light Mandalas Website:** www.light-mandalas.com
- **Official:** info@light-mandalas.com
- **Instagram:** [@mandalawave9](https://www.instagram.com/mandalawave9)
- **YouTube Channel:** Light Mandalas Official
- **Technical Support Email:** cservice@light-mandalas.com

Local agents contacts:

- **Japan:** japan@light-mandalas.com
- **Taiwan:** lovelight777888@gmail.com
- **Indonesia:** lightseed2020bali@gmail.com
- **West Malaysia:** lightworker@healingtech.org
- **East Malaysia:** malaysia@light-mandalas.com
- **Hungary:** horesnyibela@light-mandalas.com
- **Switzerland, Germany, Austria, and the Netherlands:** info@light-mandalas.ch
- **France, Spain:** france-spain@light-mandalas.com
- **USA:** usa@light-mandalas.com
- **New Zealand, Australia:** tachyonhealing.nz@gmail.com